



WorldPoint R-CAT Window for STEMI 3.25? x 6.75?

Description

- Illustrations of common EKG changes associated with Myocardial Ischemia, Injury and Infarction using a novel color-coded format
- Unique calibrated tools for evaluating heart rate, intervals (PR, QRS and QT) and significant Q-waves
- Can be used for both 12 lead EKG analysis and EKG rhythm strips
- Laminated coating which can be wiped clean with a disinfectant
- Replaces R-CAT01
- A flexible plastic pocket card with a see-through viewing window for easily comparing the patient's EKG to an established baseline.
- Illustrations of common EKG changes associated with Myocardial Ischemia, Injury and Infarction using a novel color-coded format
- Unique calibrated tools for evaluating heart rate, intervals (PR, QRS and QT) and significant Q-waves
- Can be used for both 12 lead EKG analysis and EKG rhythm strips
- Laminated coating which can be wiped clean with a disinfectant

Date

17 Jun 2025

Date Created

18 set 2022